An inter lab comparison of cyclic siloxanes in codfish collected from the Oslo Fjord

Jeremy Durham¹, Henriette Leknes², Darren Huff¹, Reinhard Gerhards³, Thomas Boehmer³, Martin Schlabach², Norman Green⁴, Roy Campbell¹, and Dave Powell¹.

¹Dow Corning Corporation, Midland, MI, USA ²Norwegian Institute for Air Research, Kjeller, Norway. ³Evonik Goldschmidt GmbH, Essen, Germany ⁴Norwegian Institute for Water Research, Oslo, Norway

Introduction


Cyclic volatile methylsiloxanes are high-volume chemicals used in a number of industrial applications and consumer products. They are currently undergoing risk assessment by the European Commission. Detection of cVMS in the Nordic environment have been reported in sediment, biota, WWTP influent and air1-3. In order to evaluate the analytical procedures for cVMS for the first time, an analytical comparison of codfish livers was performed across three separate labs: Norwegian Institute for Air Research (NILU), Evonik, and Dow Corning (DCC). Seventeen whole Atlantic cods were collected from the inner Oslofjord by the Norwegian Institute for Water Research (NIVA) and stored frozen. Each Laboratory received five to six whole frozen fish and processed the fish according to each laboratory's protocol. Portions of the liver after dissection and homogenization were divided and sent to all 3 labs for analysis. Each lab used their standard analytical method to analyze the liver homogenates for hexamethylcyclotrisiloxane (D3), octamethylcyclotetrasiloxane (D4), decamethylcyclopentasiloxane (D5) and dodecamethylcyclohexasiloxane (D6) by GC-MS. The observed concentrations from the 3 analytical laboratories were tabulated and compared for consistency.

Sampling Site Vestfjorden sampling station **VEAS** Bekkelaget Oslo **WWTP WWTP** vesodd-Kolbotn tangen Slemme stad Langhus Width: 1 km Båtstø Depth threshold Drøbak As Svelvik Tofte Atlantic cod, Gadus morhua

Study Design

- •Seventeen whole Atlantic cod fish were collected In December 2007 from the inner Oslofjord.
- •The samples were divided into three sets using a randomized block design with six blocks: 2 size groups across 3 labs.
- •Each laboratory received five to six whole frozen fish and processed the fish (harvesting and homogenization of liver) according to each laboratory's protocol.
- •At each laboratory, liver homogenates were divided and portions were sent to the other 2 labs for analysis.
- •Each lab used their standard analytical method to extract and analyze the liver homogenates

Sample preparation and analysis



Parameter	DCC	Evonik	NILU
Instrument	Agilent 6890 GC	Agilent 6890 GC	Agilent 5890 N GC
	Agilent 5973 MSD	Agilent 5973 MSD	Waters Autospec-V Ultima HRMS, R>10 000
Inlet conditions	Splitless	Cool on-column	Splitless
	1 pt.	1 μL	1 μL
	150°C		200°C
Column	Zebron ZB-5 (30 m \times 0.25 mm \times 0.25 μ m)	Agilent DB5-HT (30 m = 0.25 mm = 0.10 μm)	Agilent J&W Ultra-2 (25 m × 0.20 mm × 0.11 μm)
Oven start temp	50 °C (hold 3 min)	40°C	35°C (hold 3 min)
Transfer line	280°C	300 °C	270 °C
Analytem/z	207 (D3), 281 (D4), 355 (D5), 429 (D6)	207 (D3), 281 (D4), 355 (D5), 341 (D6)	207.0329 (D3), 281.0517 (D4), 355.0705 (D5), 429.0893 (D6
ISTD m/z	281 (M4O), 285 (BC-D4), 360 (BC-D5), 435 (BC-D6)	285(IJC-D4), 360(IJC-D5), 345(IJC-D6)	

Calculation of Detection Limits

Each laboratory used a separate definition of Limit of Detection:

DCC and Evonik: LOD = 3*SD of Matrix blank response

NILU: LOD = 3* Hexane blank response (If no blank response: Peak-to-peak S/N 3:1 (MassLynx))

Method Detection Limit:


MDL = t * SD of matrix


- DCC non-spiked codfish liver (n = 3)
- Evonik blank wolffish liver (n = 6)
- NILU spiked codfish liver (n = 5)

Summary of cVMS Concentrations

Fish ID	Processor	D3 (ng/g ww)		D4 (ng/g ww)		D5 (ng/g ww)			D6 (ng/g ww)			
		DCC	NILU	DCC	Evonik	NILU	DCC	Evonik	NILU	DCC	Evonik	NILU
OCF-07	DCC	0.1	2.9	52	70	10	379	384	62	24	33	5.9
OCF-08	DCC	0.9	3.6	66	49	42	778	743	520	67	74	44
OCF-10	DCC	0.3	3.4	94	82	77	1477	1698	1221	119	120	81
OCF-02	DCC	2.8	5.6	128	130	114	3137	3141	2790	396	387	282
OCF-15	DCC	-0.1	4.3	120	114	104	1456	1489	1238	204	195	140
OCF-13	DCC	0.7	2.1	17	-0.2	5.4	111	116	36	-6	2.2	1.6
OCF-09	Evonik	-0.1	2.8	55	28	38	433	368	347	47	55	42
OCF-06	Evonik	1.9	5.5	95	77	71	638	591	495	149	132	99
OCF-12	Evonik	7.9	13	221	280	225	1785	1938	1865	223	242	197
OCF-14	Evonik	2.2	6.0	136	131	145	2770	3023	3007	121	122	116
OCF-11	Evonik	3.0	15	103	82	105	2788	2921	2999	174	177	184
OCF-01	NILU	n.a	3.9	n.a	9.0	7.9	n.a	252	40	n.a	3.7	5.6
OCF-03	NILU	2.1	6.0	100	63	59	707	703	462	89	83	52
OCF-05	NILU	0.6	2.6	134	131	68	2009	2199	1078	138	136	62
OCF-04	NILU	4.0	5.5	103	100	110	1212	1185	1419	77	78	88
OCF-17	NILU	n.a	11	n.a	43	62	n.a	1081	1106	n.a	29	39
OCF-16	NILU	1.0	4.0	76	68	84	1542	1664	1697	171	170	154
Method Detection Limit (MDL)		7.7	3.1	3.2	7.6	4.7	23	15	8.5	11	6.5	10

The results shown are uncensored data.






Conclusions

- •While methods for extraction and processing of liver samples were quite different, overall agreement in observed concentrations was relatively good; however there were statistical differences:
- -There was no statistical difference between NILU and Evonik for D4 and D6, but there was a statistical difference for D5 -There was no statistical difference between DCC and Evonik for D4, D5 and D6
- -There was no statistical difference between NILU and DCC for D4 but there was a statistical difference for D3, D5 and D6 -No statistical difference between labs in the preparation of the liver homogenates
- •There was good agreement with previous

Oslo Fjord cod liver results 1,3

•Methods used in calculation of detection limits were different for all 3 laboratories and a limitation of this study is that a common set of data (matrix and replicates) were not defined up front to allow for a consistent determination of LOD, and MDL, across all 3 labs. This should be considered in future inter-lab comparisons.



- 1) Nordic cCouncil of Ministers, TemaNord 2005:593, Siloxanes in the Nordic Environment, L. Kaj et al., 2005. 2) Swedish Environmental Research nstitute (IVL), Report B1643, Results from the Swedish Nathonal Screening Programme 2004. Subreport 4: Siloxanes. L.
- 3) Norwegian Pollution Control Authority (SFT), TA-2269/2007. Siloxanes in the Environment of the Inner Oslofjord. M. Schlabach et.al, 2007.